

DURACELL®
ACTIVAIR®
HEARING AID

0%
MERCURY

Environmentally
friendly energy with
extra long tab



EasyTab®



Modern hearing aids need high performance batteries, Duracell Activair with EasyTab[®]

Duracell Activair hearing aid batteries work with the tried and tested zinc-air technology. The hearing aid batteries meet all requirements of modern hearing aids that have high energy requirements for additional functions and the utilisation of wireless technology.

The unique Easy Tab[®] design allows the very easy and convenient changing of Duracell Activair hearing aid batteries.



Easy to handle and easy to open.



Easy to take out and easy to hold.

Sustainably environ- mentally friendly, the contribution of Duracell Activair

Duracell Activair Mercury-free have no mercury and are environmentally friendly hearing aid batteries. These hearing aid batteries are also suitable for hearing aids that utilize wireless technology that requires extreme power.

With Duracell Activair 0% Mercury, you have a hearing aid battery characterized by:

- tried and tested technology
- long lasting performance
- constant quality and reliability
- strictest quality control tests



The unique, extra long tab
allows easy insertion.

Instructions for how to use a hearing aid battery

- Store batteries in their original packaging.
- Always store batteries at room temperature, never in the refrigerator.
- Do not remove the colored tab until you need the battery.
- Never carry loose batteries in your pocket and make sure that they never touch metallic objects.
- Always dispose of batteries in an environmentally compatible way.
- Never attempt to recharge batteries.



Your opinion is important to us

We are interested in the experiences that you have had with Duracell Activair hearing aid batteries with EasyTab[®], whether with their performance or ease of handling. Your hearing audiologist will gladly forward your opinion to us.

This overview will help you to find the right Duracell Activair battery size.

10



13



312



675



For more information, ask your
hearing audiologist.



10004.0316.2000

www.duracell-activair.com

© 2016 DURACELL, Bethel, CT 06801. Duracell and
Activair are registered trademarks of Duracell Batteries
BVBA and Duracell U.S. Operations, Inc., used under
license. All rights reserved.